


"ASB" CUSTOM FIT SYSTEM

The Adjustable Saddlebag System allows the Wander Pack to fit dogs of all sizes.




Smaller dogs:
storage pouches
close together



Larger dogs:
storage pouches
further apart

PROPER FITTING INSTRUCTIONS



1. Place sternum loop over the dogs head and rest the pack on the dogs back.
2. Adjust the stabilizer straps (located under each bag) to lower both packs so the lower part of each bag lines up with the center of the dogs side. *See illustration above*
3. Loop the two side straps under the dog's stomach and clip into the side release buckles.
4. Clip the sternum loop buckle into the clip located on the side strap.
5. Adjust sternum loop and two side straps to fit snugly. This will ensure the load is stable on the dog.
6. Load the two saddlebags with Fido's gear making sure each bag is loaded equally to prevent shifting of the bags. If you find that the pack is slipping to one side readjust the stability straps. (step 2)